ENGLISH	2
FRANÇAIS (French)	14
NEDERLANDS (Dutch)	26
ESPAÑOL (Spanish)	38
PORTUGUÊS (Portuguese)	50
DEUTSCH (German)	62
ITALIANO (Italian)	74
TÜRKÇE (Turkish)	86
JĘZYK POLSKI (Polish)	98
MAGYAR (Hungarian)	110
ΕΛΛΗΝΙΚΑ (Greek)	122
SVENSKA (Swedish)	134
ČEŠTINA (Czech)	146
БЪЛГАРСКИ (Bulgarian)	158
DANSK (Danish)	171
NORSK (Norwegian)	183
LIETUVIŲ KALBA (Lithuanian)	195
SLOVENČINA (Slovak)	207

INSTRUCTIONS FOR USE ADNEX2

ADNEX model V1.1.1

ENGLISH

Label

Device Name

Contents

ADNEX2
(01) 05430004729037 (8012) v1.1.1
Gynaia Belgium Naamsevest 88 B-3000 Leuven Belgium
Instructions for use of ADNEX2 release v1.1.1
Hyperlink to the electronic instructions for use

Symbols

Explanation of symbols

    	Unique Device Identifier = UDI-DI +UDI PI
	UDI-PI = Software version
	Date of manufacture / release date
	Manufacturer
	Consult instructions for use
	Medical Device

Device description

<i>Intended Purpose</i>	ADNEX2 is intended for use by trained healthcare professionals as an aid in the assessment and triage of women with adnexal masses. Based on selected clinical and ultrasound parameters, it provides probability estimates for outcome categories (benign, borderline, stage I invasive, stage II–IV invasive, and secondary metastatic). The outputs are intended to support, but not replace, established clinical management protocols and professional judgment. The device does not provide a diagnosis or treatment recommendation.
<i>Intended user</i>	Trained health care professionals holding current Gynaia “IOTA ovarian classification training and certification” certificate
<i>Patient target group</i>	Non-pregnant women aged 18 years and older presenting with an adnexal mass for which the Benign Descriptors are not applicable (see precautions)
<i>Contra indications</i>	Pregnancy, age < 18 year
<i>Use in combination with other devices and software</i>	N/A
<i>Limitations of use</i>	ADNEX2 is designed to support clinical decision for patient triage and is not intended for diagnostic purposes.
<i>Warnings, precautions and/or measures to be taken</i>	Use of the ADNEX2 application is contingent on the user receiving appropriate training and holding current Gynaia IOTA ovarian mass classification certification. Access to the app will not be permitted to untrained users. ADNEX2 can be used in all adnexal masses, but is intended for being used in a two-step strategy, in those cases where the Benign Descriptors do not apply. When evaluating an adnexal mass, and in the rare case you don’t have access to the ADNEX2 application due to technical issues, patient management should not be delayed, and we advise that an alternative approach be used. We recommend using IOTA Simple Rules with expert evaluation of adnexal masses when the rules cannot be applied. Information on IOTA Simple Rules and how to use them can be found on the Gynaia learning management system and in the open access publication by Timmerman <i>et al.</i> (https://doi.org/10.1002/uog.5365).

Operating instructions

Minimum System Requirements

Supported Browsers and System Requirements

Our medical device is accessed through a web application (the platform) that is compatible with a range of modern browsers on both desktop and mobile devices. This web application can also be installed on a personal device for a native-app-like experience. To ensure proper functionality, security, and the best user experience, it is essential to use a supported browser with the minimum required version. The platform will check to ensure it only works on supported devices.

Supported Browsers

The following browsers are supported for use with our medical device. Older versions of these browsers are not supported, and you must update to the specified version or higher.

- **Chrome (Desktop & Mobile):** Version 82 or higher.
- **Edge (Desktop & Mobile):** Version 82 or higher.
- **Firefox (Desktop & Mobile):** Version 77 or higher.
- **Safari (Desktop & Mobile):** Version 12 or higher.
- **Chromium-based browsers:** Any browser built on the Chromium engine (e.g., Brave, Opera) that is based on Chromium version 82 or higher should be supported, though official testing is limited to Chrome and Edge.

Browsers that have reached their "end-of-life" and are no longer receiving security updates, such as Internet Explorer, are not supported.

Hardware Requirements

The performance of the application is dependent on the hardware of your device. The following are general minimum system requirements to run the supported browsers and our application correctly on a desktop/laptop or a mobile phone.

Desktop/Laptop:

- **Processor:** A processor with support for SSE3 instructions (e.g., Intel Pentium 4 or newer) is a general requirement for modern browsers.
- **RAM:** 2GB for 64-bit browsers is the minimum, but 4GB or more is recommended for optimal performance.
- **Operating System:** An actively supported operating system such as Windows 10 or higher, macOS 12 (Monterey) or higher, or a recent version of a Linux distribution.
- **Screen Resolution:** A minimum resolution of 1024 x 768 is required.

Mobile Devices:

- **Operating System:** A modern mobile operating system that supports the latest versions of the browsers listed above (e.g., Android 10 or later, iOS 12 or later).
- **RAM:** 3GB RAM is recommended for optimal performance.

Installation

The application can be installed on a mobile device as a PWA (progressive web app) for quick access from the home screen, similar to a native app. The installation process varies slightly between operating systems.

On iPhone (iOS) using Safari:

1. Open the web application in the Safari browser.
2. Tap the **Share** button (a square with an arrow pointing upwards) at the bottom of the screen.
3. Scroll down and select **Add to Home Screen**.
4. You will be prompted to customize the name of the application. Tap **Add** in the top right corner to confirm.
5. The application icon will now appear on your home screen.

On Android using Chrome:

1. Open the web application in the Chrome browser.
2. An "Install app" prompt may appear automatically. If it does, follow the on-screen instructions.
3. If no prompt appears, tap the **three-dot menu** in the top-right corner of the browser.
4. Select **Add to Home screen** from the menu.
5. A pop-up will appear where you can customize the app's name. Tap **Add** to confirm.
6. The application icon will now be on your home screen and in your app drawer.

Performance characteristics

Clinical Benefits

ADNEX2 Clinical Benefit	Clinical Outcome Parameter
1. Improved Triage - selection for conservative management: if ADNEX2 shows a low predicted risk of malignancy, depending on the clinical context, this allows downstream guidelines to inform to follow up patients conservatively with ultrasound monitoring of the lesions	Clinical and ultrasound follow-up of patients with benign-looking lesions is safe, with a low risk of complications such as torsion, rupture or malignancy. Conservative management instead of surgery allows avoidance of surgical complications and impact on fertility
2. Improved Triage - referral to oncology centers: an increased risk of malignancy predicted by ADNEX2 will allow downstream guidelines to inform referral of the patient to oncology centers	Survival benefit
3. Improved Triage - subclassification of malignant tumors will allow downstream guidelines to inform patient management	Patients with malignant ovarian tumors will need different treatment depending on the subtype of tumor. Borderline tumors can be treated with less aggressive techniques than invasive tumors, which is of interest when fertility preservation is desired. Furthermore, stage I ovarian cancer may be managed more conservatively than stage II-IV disease. Finally, for cancers of other primary origin metastasized to the ovary, treatment depends on the primary tumor

Clinical Performance

ADNEX2 Clinical Performance	Performance endpoint/measure
1. ADNEX2 is able to give probabilities to distinguish between benign and malignant masses	The AUC for discrimination between benign and malignant masses of ADNEX2 without CA125 and with CA125 if available is 0.95 (0.94 to 0.96). The same AUC is obtained when used as part of the two-step strategy.
2. ADNEX2 is able to give probabilities to distinguish different subclasses of malignant tumors (borderline, stage 1 invasive, stage 2-4 invasive and secondary metastatic tumors)	AUCs between malignant subtypes varied between 0.80 and 0.98 for ADNEX2 with CA125 if available. For ADNEX2 without CA125, AUCs between malignant subtypes varied between 0.65 and 0.97.
3. The risks predicted by ADNEX2 both with and without CA125 are well calibrated ¹	ADNEX2 showed close to perfect calibration in the whole population but calibration in individual centers varies.

¹ **Calibration:** the degree to which the risk probabilities predicted by ADNEX2 correspond to the actual observed outcomes in patients.

Maintenance

No maintenance is required.

For both security and technical performance reasons (including speed, user experience, and accurate risk visualization) it is strongly recommended to use the latest available version of the operating system and browser on the device running the ADNEX2 application.

Risks - Information for Safety

No risks are associated with installing and operating the software.

	Incorrect data entry, including misinterpretation of ultrasound features or incorrect entry of clinical variables may potentially impact the risk estimation
	Use by untrained users, potentially leading to incorrect entry of clinical variables and/or misinterpretation of results
	Inappropriate use in populations where the application has not been validated, namely patients under the age of 18 and pregnant patients. Risks may also increase if users bypass or ignore warnings and applicability checks embedded in the software
	Cybersecurity risks, such as unauthorized access or misrepresentation of user identity (e.g., accessing the software using another person's credentials). These may be mitigated by access controls and using secure, up-to-date browsers and devices.
	The ADNEX2 risk calculation is only indicative to support physicians in selecting patients for appropriate management within established management protocols. It is not intended to be used to establish a diagnosis.
	Even though ADNEX2 has a high level of performance, users should be aware that false positive and false negative results can rarely occur
	Calibration performance is very heterogenous and is based on a large multicenter study. Performance in individual centers can differ from the reported calibration.

Data Security

Area	Physician Responsibility (Duty of Care)
Active Use & Logout	NEVER leave the application logged in and unattended. After every consultation or period of use, you must securely log out (or rely on the system's <i>immediate</i> automatic log-off feature).
Screen Privacy (Public View)	ENSURE Visual Privacy while using the application. Position your screen (PC, smartphone, tablet) so that it cannot be viewed by unauthorized persons.
Device Security	Your device must be secured with a strong PIN, password, or biometrics (fingerprint/face ID). Ensure the device's operating system and the application are updated to the latest security versions.
Data Transmission	Only use the application over trusted, encrypted networks (e.g., your clinic's secured Wi-Fi or a trusted VPN).
Personal Data Export	Do not download or take screenshots of Protected Health Information (PHI) to your personal device's local storage unless absolutely necessary and immediately secure or delete the file after use.

DISCLAIMER

The manufacturer is not responsible for security breaches, data loss, or system compromise resulting from the user's failure to adhere to the security and data protection responsibilities outlined in this IFU, including but not limited to, the use of weak passwords, the use of unapproved operating environments, or failure to apply recommended software updates.

Instructions for Use for this product are provided electronically (eIFU). No paper Instructions for Use are supplied by default.

Paper copies of the Instructions for Use are available upon request. Printed copies are considered uncontrolled documents. Users are responsible for verifying that the version in use is the current and applicable version.

Manufacturer

Contact details:

Gynaia Belgium
Naamsevest 88
B-3000 Leuven
Belgium
T: +32 497 79 89 99
E: support@gynaia.com
SRN: BE-MF-000048109

United Kingdom – UK Responsible Person

Qarad UK Ltd
8 Northumberland Ave
Westminster
London WC2N 5BY
United Kingdom

ukrp@qbdgroup.com

Registration number: 13004437
VAT Number: GB364 2759 75

Incidents and Support

For technical support or questions contact Gynaia support team via support@gynaia.com.

Serious incidents should be reported immediately to Gynaia and to the competent authority of the Member State in which the user/patient is established or the incident took place.

Document Version

Version	Issue date	Description of modification
1.0	27 October 2025	New document
1.1	28 November 2025	Updated version information because of a software patch to expose ADNEX base relative risk scores.
2.0	18 December 2025	Added clinical Benefits and NB number to CE label. Added Data security section and disclaimer.

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